


SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: ALKYL GLYCOSIDES

HPLC



- 1: Methyl α , β -D Glucopyranoside
- 2: Hexyl β -D Glucopyranoside
- 3: Octyl β -D Galactopyranoside
- 4: Octyl β -D Glucopyranoside
- 5: Octyl β -D Xylopyranoside
- 6: Decyl β -D Maltoside
- 7: Decyl β -D Glucopyranoside

CONDITIONS

| | |
|-----------------|--|
| Column | μ Bondapak Phenyl -10 μ m 30 x 4.6 mm |
| Injected Volume | 10 μ l (Injected Quantity 4 μ g) |
| Eluent | Acetonitrile/Water (20/80) |
| Flow Rate | 1.5 ml/min |
| Gas | Air |
| Pressure | 2.5 Bars |
| Temperature | 40° C |

Analyzed By : LCBA 45100 Orleans FRANCE