


SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: ALKYL GLYCOSIDES

SFC

- 1: Octyl β -D Galactopyranoside
- 2: Octyl β -D Glucopyranoside
- 3: Octyl β -D Xylopyranoside



CONDITIONS

Column	<i>Zorbax Phenyl - 7 μm 15 x 4.6 mm</i>
Eluent	<i>CO₂/Methanol (95.2/4.8)</i>
Flow Rate	<i>1.5 ml/min</i>
Gas	<i>Air</i>
Pressure	<i>159 Bars</i>
Temperature	<i>Detector 40° C, Column 40° C</i>

Analyzed By : LCBA 45100 Orleans FRANCE