


SEDEX Low Temperature Evaporative Light Scattering Detector


SAMPLE: MALTODEXTRINS

HPLC
Gradient

SEDEX Evaporative Light Scattering Detector can detect concentrations of 1000 ppm of DP 1 to DP 17



Refractive Index detector can detect concentrations of 50,000 of DP 1 to DP 10



Gradient: 60% to 80% Acetonitrile in 20 min

CONDITIONS

Column	<i>Erbasil NH₂</i> <i>25 x 4.6 mm</i>
Injected Volume	<i>20 µl, Injected Quantity 20 µg</i>
Eluent	<i>A: Acetonitrile, B: Water</i>
Flow Rate	<i>1.5 ml/min</i>
Pressure	<i>2.5 Bars</i>
Gas	<i>Air</i>
Temperature	<i>30° C</i>