


SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: ALKYL GLYCOSIDES

HPLC



- 1: Methyl α , β -D Glucopyranoside
- 2: Hexyl β -D Glucopyranoside
- 3: Octyl β -D Galactopyranoside
- 4: Octyl β -D Glucopyranoside
- 5: Octyl β -D Xylopyranoside
- 6: Decyl β -D Maltoside
- 7: Decyl β -D Glucopyranoside

CONDITIONS

Column	μ Bondapak Phenyl -10 μ m 30 x 4.6 mm
Injected Volume	10 μ l (Injected Quantity 4 μ g)
Eluent	Acetonitrile/Water (20/80)
Flow Rate	1.5 ml/min
Gas	Air
Pressure	2.5 Bars
Temperature	40° C

Analyzed By : LCBA 45100 Orleans FRANCE