


Assessment of a New Generation of Evaporative Light-Scattering Detectors for Liquid Chromatography: Sensitivity, Linearity, Dynamic Range, Analyte Dispersion and Response Variation with Eluent Composition



Eric VERETTE, Ph.D.
SEDERE S.A.S., France



Abstract:


- Among the detectors available in Liquid Chromatography, Evaporative Light-Scattering Detector (ELSD) became in recent years a well established instrument thanks to several theoretical studies based on fundamental investigations and numerous applications provided during the last thirty years. Indeed, ELSD is considered as a nearly Universal, powerful and cost-effective technique, and is ideally appropriate for a great majority of liquid chromatography applications containing chromophoric and non-chromophoric compounds. However, some limitations have been attributed to this type of detection in the past, such as low sensitivities, non-linearity, reduced dynamic range, large analyte dispersion and wide response variation with the mobile phase composition.

- Today, a new ELSD model is proposed which offers a genuine and efficient Low-Temperature technology (LT-ELSD™) combined with an innovative detection chamber. The overall design of this new detector results in a significant increase of sensitivity providing limits of detection down to the sub-nanogram levels even for semi-volatile compounds. Moreover, it provides an improved direct linearity with correlation coefficients reaching 0.999 and an extended overall dynamic range exceeding the four orders of magnitude. Also, this model is optimized for the more recent U-HPLC technique giving peak widths below 1 second. Finally, a study on the response variation with the mobile phase composition showed a smaller effect compared to other aerosol-based detectors.

- In this work, several examples and data are provided on these topics to emphasize the great advancement in this technology, thus bringing to the analyst a relevant and cost-effective solution to their chromatography challenges.

I - Sensitivity


Application: Generic HPLC/LT-ELSD Method for Lipids



This example shows very High Sensitivities obtained with a real HPLC/LT-ELSD application. LODs are much below 10ng on column for all compounds (except for Lauric acid which has got a high semi-volatility feature), and even at the Picogram Levels for other semi-volatile compounds such as Fatty alcohols and some Fatty acids.

II - Linearity


Application: HPLC/LT-ELSD Analysis of Phenolic Acids



This application shows very good Direct Linearities for the 5 Phenolic acids with concentrations from 50ng to 10000ng on column. Correlation coefficients ranged from 0.9986 to 0.9999. Additionally, repeatability in this application varied from 0.6% to 2.5%. Elsewhere, the figure also presents a small variation of the compound responses with identical concentrations and gave a relative variation percentage value of +/- 8% (e.g.: Gentisic acid and Protocatechuic acid curves are overlapping). This demonstrates that SEDEX 90LT, as an efficient and genuine mass detector, provides uniform responses.


III - Dynamic Range

Application: HPLC/LT-ELSD Analysis of Hydrocortisone



IV - Efficiency, Peak Widths

Application: U-HPLC/LT-ELSD Analysis of Phenolic Acids




Courtesy of D. Guilleme, University of Geneva, Switzerland.

These 2 applications demonstrate the good suitability of SEDEX 90LT for U-HPLC. In the application on Phenolic acids, peak widths (1/2 height) ranged from 0.8 to 1.0 second.

V - Response Variation with H2O/ACN Composition

Application: HPLC/LT-ELSD Analysis of Hydrocortisone



Conclusion:

The applications developed above clearly show the advantages of the new SEDEX 90LT ELSD and particularly in regards to:

- Sensitivity with low nanogram and even sub-nanogram levels,
- Linearity with correlation coefficients reaching 0.999,
- Dynamic range of 4 orders of magnitude,
- Excellent suitability for U-HPLC with peak widths below 1 second,
- Low response variation compared to all other aerosol-based detectors.

This work also demonstrates the significant advancements of the new SEDEX Evaporative Light-Scattering Detector resulting from the combination of an efficient and genuine Low Temperature technology and an innovative detection chamber. These outstanding new features offer now to the analyst a Universal, powerful, versatile and cost-effective solution to their separation and quantification challenges.