



**Sensitivity  
Flexibility  
Experience**

# A New Generation of Evaporative Light-Scattering Detectors for Liquid Chromatography: Universality, High Performance and Robustness in Pharmaceutical Analysis - An Application Review in HPLC and U-HPLC

**Eric VERETTE, Ph.D.**  
**SEDERE S.A.S., France**

## Abstract


Among the detectors available in Liquid Chromatography (LC), Evaporative Light-Scattering Detector (ELSD) became in recent years a well established instrument thanks to several theoretical studies based on fundamental investigations and numerous applications provided during the last thirty years. Indeed, ELSD is considered as a nearly Universal, powerful, reliable and cost-effective technique, and is ideally appropriate in Pharmaceutical industry for a great variety of LC applications containing chromophoric and non-chromophoric compounds.

Today, an ELSD model based on a recent and unique concept is proposed which offers a genuine and efficient Low-Temperature technology (LT-ELSD™) combined to an innovative detection chamber. The overall design of this ultimate detector results in a significant increase of sensitivity providing typical limits of detection down to the very low nanogram levels for non-volatile and semi-volatile compounds. It provides an improved overall direct linearity with correlation coefficients over 0.99, consistent responses independent of the analytes chemical structure and an extended dynamic range exceeding the four orders of magnitude (from low ng to high µg levels on column). Also, this model is optimized for the recent U-HPLC technique giving peak widths of less than 1 second.

To show the strength and the versatility of this ELSD model, several relevant LC applications in Pharmaceutical analysis are developed in this work. These applications use the most recent LC media, such as multi-mode, HILIC and sub-two-micron or fused-core particle phases, allowing outstanding separations and simultaneous analyses of a wide range of chemical and biochemical compounds.

The topics presented here are focused on:


- Sensitivity and robustness in LC on the full pH range.
- Rapid and simultaneous separation of API, counterions and impurities.
- Response consistency compared to UV detectors at several wavelengths.
- Relevance in the determination of both chromophoric and non-chromophoric solutes in the analysis of natural products or TCM, such as Ginkgo Biloba.
- Simplified and cost-effective alternative in the analysis of non-chromophoric compounds such as aminoglycoside antibiotics, thus avoiding any tedious derivatization step, ion pair reagents and specific detectors.



**HPLC / LT-ELSD System**

## I - Sensitivity


Application: Generic HPLC/LT-ELSD Analysis of Several Lipid Compounds



The results show very high sensitivities obtained in a real HPLC/LT-ELSD application. Obtained LODs are much below 10ng on column for all compounds (except for Lauric acid which is characterized by a high vapor pressure), and even at the Picogram Levels for some other semi-volatile compounds belonging to the groups of fatty alcohols and fatty acids.

## II - API and their Counterions


Application: Global HPLC/LT-ELSD Method for the Simultaneous Analysis of Polar and Non-Polar, Neutral, Acidic and Basic Pharmaceutical Drugs and their Respective Counterions



This chromatogram shows that the combination of an efficient multimodal stationary phase and a single Universal ELS detector allows the quick and easy simultaneous determination of a wide range of compounds characterized by different polarities, including inorganic anions and cations.

## III - Impurity Assessment

Application: Simultaneous HPLC/LT-ELSD Analysis of Imipramine, its Counterion and an Impurity



In this example SEDEX 90LT detects an impurity (sodium) at a level of 0.05% of the major compound (Imipramine). With this detector, there are 12 gains available, which means that this percentage could typically go down to less than 0.01% (more than 4 orders of magnitude for the overall dynamic range).

## IV - Response Consistency

Application: Fast HPLC/LT-ELSD/DAD Analysis of Non-Volatile Compounds with Different Chemical Structures

