



SEDEX Low Temperature Evaporative Light Scattering Detector

SAMPLE: PREDNISOLONE (with Impurities)

SFC



CONDITIONS

Column	<i>Zorbax TMS - 7 μm 250 x 4.6 mm</i>
Eluent	<i>CO₂/MeOH (89.3/10.7)</i>
Flow Rate	<i>4.6 ml/min</i>
Pressure	<i>3000 psi</i>
Temperature	<i>Detector: 40° C, Column : 40° C</i>

Analyzed By : LCBA 45100 Orleans FRANCE

S.E.D.E.R.E.
PARC VOLTA - BP 27
9, rue Parmentier
94141 Alfortville Cedex FRANCE

Application #13